

No. 91578UT



Oxysept 1 Step® Neutralizing Tablet contains HPMC (Hypermellose), a built-in artificial tear moisturizing lubricant which gives the neutralized solution characteristics of natural tears. HPMC forms a hydrophilic lubricating layer surrounding your lens to keep your lens moist for unsurpassed comfort.

The Outer HPMC layer of Oxysept 1 Step® Neutralizing Tablet enables delayed neutralization, effecting the neutralization action after the solution has completely killed all harmful bacteria and microbes (even Acanthamoeba) in 6 hours.

Oxysept 1 Step® Neutralizing Tablets contains Vitamin B12, an unique color indicator. The Oxysept 1 Step® Neutralizing Tablet turns the disinfection solution light rosy pink as it dissolves, thus ensuring 100% neutralizing in 6 hours. By using a new tablet everyday, there is no need to worry about the reduced effectiveness of neutralization unlike using an ageing neutralizing disc. Also, Vitamin B12 is a naturally occurring vitamin in our body.



Manufactured by  
Nicholas Piramal India Limited  
(889, Plot 67-70, Sector II, Pitampura 454 775  
Madhya Pradesh, India)  
Manufactured for Advanced Medical Optics  
(2611, 1700 E. St. Andrew Place, Santa Ana  
California 92709 USA)

XXXXXXXXXX



4 897004 110073

Oxysept 1 Step® Neutralizing Tablets

AMO  
**Oxysept**  
**1Step®**

Neutralizing Tablets



**Contains Vitamin B12  
colour indicator**

**For all soft contact lenses**

**12 Tablets**

Oxysept 1 Step® Neutralizing Tablets